



DC Bead™

Ordering Information

Label Colour	Nominal Bead Size	Volume of Beads	Product Code
Yellow	100 – 300 µm	2ml	DC2V103
Blue	300 – 500 µm	2ml	DC2V305
Red	500 – 700 µm	2ml	DC2V507
Green	700 – 900 µm	2ml	DC2V709

10ml glass vial. Each vial contains approximately 2ml of DC Bead™ in non-pyrogenic, sterile, physiological buffered saline. Total volume of saline and DC Bead™ is approximately 8ml. The vial is stopper sealed by an aluminium cap equipped with a colour-coded lid.

PRECISION TACE
Less More Every Time